HEAD MEASUREMENT INSTRUCTIONS

STEP 1 - LOCATE THE MID POINT

Run an imaginary line across from behind one ear all away up and across the dogs head to the back of the other ear you will find the mid-point located in the centre of this line. This mid-point is used to establish the HEAD MEASUREMENT.

STEP 2 - ESTABLISH THE HEAD MEASUREMENT

Measure from the mid-point to the tip of the dogs nose

- DO NOT follow the contours of the dogs face.
- Measure in a **STRAIGHT LINE ONLY** (you may need to lift a straight ruler away from the dogs head see below diagrams)

The resulting measurement represents the HEAD MEASUREMENT you will use to size your OPTIVIZOR using the table provided.

FINAL SIZING

Using the tables provided size the OPTIVIZOR by HEAD MEASUREMENT AND by CANINE WEIGHT. You may find you have 2 different sizes according to the 2 charts. This is because, as we know, dogs come in all shapes and sizes. Simply pick the smaller size. It is always better to size down when in doubt.

A good fit is critical function correctly the best chance of accepting this alternative to the dreaded "CONE OF SHAME"

When fitting your new OPTIVIZOR make sure the neck collar and chin strap is done up firmly otherwise the dog is less likely to accept the collar. Depending on the temperament of your dog, you may loosen these slightly once the OPTIVIZOR has been accepted by

SIZE MATTERS!

for the OPTIVIZOR to and to give your dog

A SNUG FIT!

your special friend.

4.3 - 5 5 - 6 6-7 7-8.3 8.3-9.4 9.5+ inches (Equivelent E-Collar size in Brackets) TOY (10) (P) **SMALL (12.5) SMALL MEDIUM (12.5 - 20) MEDIUM (20 - 25)** LARGE (25 - 30) X LARGE (30 - 40)

PROVIZOR

OPTIVIZOR DIAGRAMS





